

# Transport Chairs

## Features and Benefits

### Aluminum Transport Chairs

- Aluminum construction with a 19-inch seat is lightweight and strong, making it easy to transport patients up to 300 lbs
- Black, nylon upholstery is breathable, comfortable, durable and easy-to-clean
- Full, fixed arms are padded for patient comfort
- Tool-free, swing-away footrests remove by releasing the lock, swinging outward and lifting off to accommodate user preferences
- Foldable back assists in transporting patient to and from seated position
- Seat belt buckles for added transport safety
- Solid, 8-inch caster wheels with rear locks are maintenance-free and provide a smooth, stable ride over most surfaces and security in parked position
- Only 19 lbs and easily collapses in two-step method for convenient storage and travel
- Plastic, ergonomic handles for attendant's comfort and grip
- Not made with natural rubber latex



## ORDERING INFORMATION

TRANSPORT CHAIRS				
MFR #	DESCRIPTION	WEIGHT CAPACITY	HCPCS	QTY
146-TR39E-SV ⚠	Lightweight Silver Vein Steel	250 lbs	E1038	1 ea
146-ATC19-BL ⚠	Lightweight Blue Aluminum	300 lbs	E1038	1 ea
146-RTLFW19RW-RD ⚠	Lightweight Red Aluminum	300 lbs	E1038	1 ea

**⚠ WARNING:** This product can expose you to chemicals known to the State of California to cause birth defects or other reproductive harm. For more information, go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).