

Rollators

Features and benefits

Durable rollators

- Steel frame construction safely supports up to 300 lbs
- Designed with folding hinge for easy storage and portability
- Ergonomically-designed handles are easy-to-grip, reduce hand pressure and are adjustable to accommodate users' height and preferences
- Lever allows for independent handle height and angle adjustment
- Oversized padded seat and backrest are made of strong polyurethane material for durability, have a contoured design for ergonomic comfort, and provide cut and split resistance to prevent rubbing and irritation
- Hinged padded backrest can be folded up, down or removed
- Deluxe loop-lock system operates like a bicycle hand brake to slow down or temporarily stop and easily locks into stationary position by pushing handle into down-position
- Swiveling, 7 1/2-inch, non-marking casters are ideal for indoor and outdoor use
- Includes convenient cloth storage pouch for personal items
- 12-inch seat width; 21-inch seat-to-floor height
- Not made with natural rubber latex
- HCPCS: E0143, E0156

