

## Transport Chairs

When assisted by a caregiver, transport chairs merge the functionality of a wheelchair with the portability of a rollator. McKesson transport chairs are lightweight, easy to use, comfortable and collapsible for convenient storage. Available in steel or aluminum construction, McKesson transport chairs feature built-in seat buckles and rear wheel locks for added safety.

### Features and Benefits

#### Steel Transport Chairs

- Durable steel construction with an attractive silver vein finish adds scratch-resistance and supports up to 250 lbs
- Black, nylon upholstery is breathable, comfortable, durable and easy-to-clean
- Full, fixed arms are padded for patient comfort
- Tool-free, swing-away footrests remove by releasing the lock, swinging outward and lifting off to accommodate user preferences
- Foldable back assists in transporting patient to and from seated position
- Seat belt buckles to enhance passenger safety
- Solid, 8-inch caster wheels with rear locks are maintenance-free and provide a smooth, stable ride over most surfaces and security in parked position
- Features folding back and upholstery handles to conveniently collapse in simple two-step method to easily store or transport the lightweight 26 lb chair
- Plastic, ergonomic handles for attendant's comfort and grip
- Not made with natural rubber latex



#### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

#### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.